---

# EIAgen

## Detect HIV 4 Total Screening Kit

---

**REF** 081311  96

**REF** 081312  192

**REF** 081315  480

**IVD**

This package insert must be read carefully before product use.  
Package insert instructions must be carefully followed.

Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this package insert.

**Manufacturer:**  
**Adaltis S.r.l**  
Via Durini, 27  
20122 Milano (Italy)  
Tel. +39-0774-5791 - Fax +39-0774-353085  
[www.adaltis.net](http://www.adaltis.net)

**en**

**SYMBOLS USED ON LABELS**

<b>English</b> <b>EN</b>							
	In Vitro Diagnostic Medical Device	Catalogue Number	Lot Number	Attention, See Instructions For Use	Temperature Limitation	Use By	Number of Test
					<b>MICROPLATE</b>	<b>DIL SPE</b>	<b>CONTROL -</b>
	Manufacturer	Keep away from Sunlight	Biological Risk	Date of Manufacture	Microplate	Sample Diluent	Negative Control
	<b>CONTROL + HIV-1</b>	<b>CONTROL + HIV-2</b>	<b>CONTROL + HIV p-24</b>	<b>CONJ 1</b>			
	Positive Control HIV-1	Positive Control HIV-2	Positive Control HIV p-24	Conjugate # 1			
	<b>CONJ 2 100X</b>	<b>CONJ 2 DIL</b>	<b>SUBS TMB</b>	<b>SOLN STOP</b>	<b>WASH BUF 25X</b>		
	Conjugate # 2 Concentrate 100x	Conjugate # 2 Diluent	Substrate TMB	Stop Solution	Wash Buffer Concentrate 25x		

**A. INTENDED USE**

The EIAgen Detect HIV 4 Total Screening Kit is a 4th generation solid phase Enzyme-Linked Immunosorbent assay using a mixture of antigens and antibodies for the *in vitro* diagnostic screening in human serum or plasma (EDTA, Heparin and Citrate) of antibodies to HIV-1, HIV-2 and HIV-1 p24 antigen.

This kit is a combined Ag/Ab assay and is not to be used for the detection of HIV-1 p24 antigen alone. This kit is for *in vitro* diagnostic use by a health-care professional and will not be sold to the general public.

**B. INTRODUCTION**

Acquired immune deficiency syndrome (AIDS) is a set of symptoms resulting from the incapacitation of the human immune system caused by the Human Immunodeficiency Virus (HIV). HIV infection may progress to a symptomatic phase that is characterized by opportunistic infections and that may cause death.

The etiological agent of AIDS, HIV, targets specific types of T cells causing Lymphopenia and affecting T cell mediated immunity. HIV is a member of a retrovirus family with two sub-families: HIV-1 and HIV-2. HIV-1 is more virulent and transmittable than HIV-2. HIV-1 is the cause of HIV infections globally, whereas HIV-2 is found predominantly in the countries of West Africa. As serological cross-reaction between HIV-1 and HIV-2 is highly variable and dependent on the tested sample, antigens for the specific detection of both HIV-1 and HIV-2 are included in the assay.

HIV is transmitted through sexual contact with infected persons, sharing needles and syringes with infected people and transfusion of contaminated blood. Enzyme-Linked Immunosorbent Assays (such as the EIAgen Detect HIV 4 Total Screening ELISA 4th gen) are recommended for screening human blood and plasma for the presence of anti-HIV antibodies and HIV-1 p24 antigen. The presence of anti-HIV-1 and/or anti-HIV-2 in the blood indicates potential infection with HIV-1 and/or HIV-2 and consequently this blood should not be used for transfusion or for manufacture of injectable products.

**C. PRINCIPLE OF THE TEST**

Antigens representing epitopes of HIV-1 gp41 and HIV-2 gp36 are coated onto microplate wells together with monoclonal antibodies against HIV-1 p24. Serum or plasma sample is added to the well and if antibodies specific for HIV-1 and/or HIV-2 (IgG, IgM or IgA) are present in the sample, stable complexes will be formed with the HIV antigens attached to the well. HIV-1 p24 antigen, if present will bind simultaneously to the antibodies in the well and to the detector antibodies present in the Sample Diluent. Non-reactive antibodies are removed by washing. Stable antigen-antibody complexes are identified through the successive addition of biotinylated antigens and horseradish peroxidase (HRP) conjugated streptavidin. These antibody-antigen complexes are quantified through the catalytic activity of horseradish peroxidase. Peroxidase substrate solution is added and is converted to a blue-coloured product. A positive sample generates a dark blue colour while faint blue colour or colourless wells indicate a negative sample. Upon adding stopping solution, the colour of the

solution will change from blue to yellow. Optical Density (OD) is measured with a spectrophotometer (ELISA reader) at 450nm with reference wavelength at 600-650nm and is in proportion to the amount of anti-HIV1/2-antibodies and HIV-1 p24 present in the sample.

**D. COMPONENTS**

The kit contains reagents for 96 tests (code 081311), or 192 tests (code 081312) or 480 tests (code 081315).

<b>Microplate</b>	1
<b>Negative Control</b>	1x2 mL/vial
<b>Positive Control HIV-1</b>	1x2 mL/vial
<b>Positive Control HIV-2</b>	1x2 mL/vial
<b>Positive Control HIV p-24</b>	1x2 mL/vial
<b>Conjugate # 1</b>	1x25 mL/vial
<b>Conjugate # 2 Concentrate 100x</b>	1x0.25 mL/vial
<b>Conjugate # 2 Diluent</b>	1x25 mL/vial
<b>Sample Diluent</b>	1x12.5 mL/vial
<b>Substrate TMB</b>	1x40 mL/vial
<b>Stop Solution</b>	1x15 mL/vial
<b>Wash Buffer Concentrate 25x</b>	1x50 mL/vial
<b>Plate sealing foils</b>	2
<b>Number of tests</b>	96
<b>Code</b>	081311

<b>Microplate</b>	2
<b>Negative Control</b>	1x2 mL/vial
<b>Positive Control HIV-1</b>	1x2 mL/vial
<b>Positive Control HIV-2</b>	1x2 mL/vial
<b>Positive Control HIV p-24</b>	1x2 mL/vial
<b>Conjugate # 1</b>	2x25 mL/vial
<b>Conjugate # 2 Concentrate 100x</b>	2x0.25 mL/vial
<b>Conjugate # 2 Diluent</b>	2x25 mL/vial
<b>Sample Diluent</b>	2x12.5 mL/vial
<b>Substrate TMB</b>	2x40 mL/vial
<b>Stop Solution</b>	1x40 mL/vial
<b>Wash Buffer Concentrate 25x</b>	2x50 mL/vial
<b>Plate sealing foils</b>	4
<b>Number of tests</b>	192
<b>Code</b>	081312

<b>Microplate</b>	5
<b>Negative Control</b>	1x4 mL/vial
<b>Positive Control HIV-1</b>	1x4 mL/vial
<b>Positive Control HIV-2</b>	1x4 mL/vial
<b>Positive Control HIV p-24</b>	1x4 mL/vial
<b>Conjugate # 1</b>	5x25 mL/vial
<b>Conjugate # 2 Concentrate 100x</b>	5x0.25 mL/vial
<b>Conjugate # 2 Diluent</b>	5x25 mL/vial
<b>Sample Diluent</b>	5x12.5 mL/vial
<b>Substrate TMB</b>	3x40 mL/vial
<b>Stop Solution</b>	2x40 mL/vial
<b>Wash Buffer Concentrate 25x</b>	5x50 mL/vial
<b>Plate sealing foils</b>	10
<b>Number of tests</b>	480
<b>Code</b>	081315

**1. Microplate**

12 strips of 8 breakable wells coated with HIV specific gp36 and gp41 peptides and gp41 protein with Monoclonal Antibodies specific to the HIV-1 p24 Ag. Plates are sealed into an aluminium pouch with desiccant.

Bring the microplate to room temperature (+18...24°C) before opening the bag. Unused strips have to be returned into the pouch and the pouch has to be sealed and stored back to +2...8°C, in presence of the desiccant.

### 2. Negative Control

Ready to use control. It contains human serum negative for HIV antibodies and for p24 antigen, and 0.05% Proclin 300 as preservative.

### 3. Positive Control HIV-1

Ready to use control. It contains human serum positive to Anti-HIV 1 and 0.05% Proclin 300 as preservative.

**Important Note:** *The absence of viable pathogens in the positive control can not be fully ensured, and therefore, the control should be handled as potentially biohazardous, in accordance with good laboratory practices*

### 4. Positive Control HIV-2

Ready to use control. It contains human serum positive to Anti-HIV 2 and 0.05% Proclin 300 as preservative.

**Important Note:** *The absence of viable pathogens in the positive control can not be fully ensured, and therefore, the control should be handled as potentially biohazardous, in accordance with good laboratory practices*

### 5. Positive Control HIV p-24

Ready to use control. It contains recombinant HIV-1 p24 and 0.05% Proclin 300 as preservative.

### 6. Conjugate # 1

Ready to use and blue colour coded reagent. It contains a mixture biotinylated antigens in phosphate buffer and 0.05% Proclin 300 as preservative.

### 7. Conjugate # 2 Concentrate 100x

100x solution concentrate. It contains Horseradish Peroxidase conjugated with streptavidin in specific stabilizer buffer with preservative.

Reagent to be diluted with the Conjugate # 2 Diluent.

**Important Note:** *Any unused portion of this diluted Conjugate # 2 Solution may be stored at +2...8°C for no more than 6 days.*

### 8. Conjugate # 2 Diluent

Ready to use and orange colour coded reagent.

Reagent is used for the dilution of the Conjugate # 2 Concentrate 100x. It contains phosphate buffer and 0.05% Proclin 300 as preservative.

### 9. Sample Diluent

Ready to use and red colour coded reagent.

It contains active components and 1.0% v/v aggregated mouse serum proteins, 0.1% w/v aggregated human IgG and 0.05% Proclin 300 as preservative.

### 10. Substrate TMB

Ready-to-use component. It contains 0.26 mg/mL of 3,3',5,5' Tetramethylbenzidine (TMB) and 0.01% w/v of Hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>), in citrate buffer.

Mix gently before use.

### Note:

*To be stored protected from light as sensitive to strong illumination.*

### 11. Stop Solution

Ready-to-use component.

It contains 0,3 M H<sub>2</sub>SO<sub>4</sub> solution.

Mix gently before use.

**The MSDS is available upon request of laboratory personnel.**

### 12. Wash Buffer Concentrate 25x

25x concentrated solution.

It contains 0.2% Proclin 300 as preservative. Once diluted, the wash solution (wash buffer diluted) contains phosphate buffer saline and Proclin 300 and Tween 20 as preservatives.

Once diluted, the wash solution is stable for 6 days at +2...8°C.

### E. MATERIALS REQUIRED BUT NOT PROVIDED

1. Calibrated Micropipettes of 200 µL, 50 µL and 10 µL and disposable plastic tips.
2. EIA grade water (double distilled or deionised, charcoal treated to remove oxidizing chemicals used as disinfectants).
3. Timer with 60 minute range or higher.
4. Absorbent paper tissues.
5. Calibrated ELISA microplate thermostatic incubator (dry or wet) set at +37°C (± 0.5°C tolerance).
6. Calibrated ELISA microwell reader with 450nm (reading) and if possible with 600-650nm (blanking) filters.
7. Calibrated ELISA microplate washer.
8. Vortex or similar mixing tools.

### F. WARNINGS AND PRECAUTIONS

1. The kit has to be used by skilled and properly trained technical personnel only, under the supervision of a medical doctor responsible of the laboratory. This package insert must be read carefully before product use.
2. The kit has to be used in a laboratory certified and qualified by the national authority in that field (Ministry of Health or similar entity) to carry out this type of analysis.
3. All the personnel involved in performing the assay have to wear protective laboratory clothes, talc-free gloves and glasses. The use of any sharp (needles) or cutting (blades) devices should be avoided. All the personnel involved should be trained in biosafety procedures, as recommended by the Center for Disease Control, Atlanta, U.S. and reported in the National Institute of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
4. All the personnel involved in sample handling should be vaccinated for HBV and HAV, for which vaccines are available, safe and effective.
5. The laboratory environment should be controlled so as to avoid contaminants such as dust or air-born microbial agents, when opening kit vials and microplates and when performing the test. Protect the Substrate TMB (TMB & H<sub>2</sub>O<sub>2</sub>) from strong light

and avoid vibration of the bench surface where the test is undertaken.

6. Upon receipt, store the kit at +2...8°C into a temperature controlled refrigerator or cold room.
7. Do not interchange components between different lots of the kits. It is recommended that components between two kits of the same lot should not be interchanged.
8. Check that the reagents are clear and do not contain visible heavy particles or aggregates. If not, advise the laboratory supervisor to initiate the necessary procedures for kit replacement.
9. Avoid cross-contamination between serum/plasma samples by using disposable tips and changing them after each sample.
10. Avoid cross-contamination between kit reagents by using disposable tips and changing them between the use of each one.
11. Do not use the kit after the expiration date stated on external (carton box/secondary container) and internal (vials) labels.
12. Treat all specimens as potentially infective. All human serum specimens should be handled at Biosafety Level 2, as recommended by the Center for Disease Control, Atlanta, U.S. in compliance with what reported in the Institutes of Health's publication: "Biosafety in Microbiological and Biomedical Laboratories", ed. 1984.
13. The use of disposable plastic-ware is recommended in the preparation of the washing solution or in transferring components into other containers of automated workstations, in order to avoid contamination.
14. Waste produced during the use of the kit has to be discarded in compliance with national directives and laws concerning laboratory waste of chemical and biological substances. In particular, liquid waste generated from the washing procedure, from residuals of controls and from samples has to be treated as potentially infective material and inactivated. Suggested procedures of inactivation are treatment with a 10% final concentration of household bleach for 16-18 hrs or heat inactivation by autoclave at 121°C for 20 min.
15. Accidental spills have to be adsorbed with paper tissues soaked with household bleach and then with water. Tissues should then be discarded in proper containers designated for laboratory/hospital waste.
16. The Stop Solution contains 0,3M sulphuric acid. Avoid contact with skin and eyes. In the event of contact, rinse immediately with plenty of water.
17. Do not smoke, eat, drink or apply cosmetics in areas in which specimens or kit reagents are handled
18. Other waste materials generated from the use of the kit (example: tips used for samples and controls, used microplates) should be handled as potentially infective and disposed according to national directives and laws concerning laboratory wastes.
19. Do not pipette by mouth.

#### **G. SPECIMEN: PREPARATION AND RECOMMENDATIONS**

1. Blood is drawn aseptically by venepuncture and plasma or serum is prepared using standard techniques of preparation of samples for clinical

laboratory analysis. No influence has been observed in the preparation of the sample with citrate, EDTA and heparin.

2. Avoid any addition of preservatives to samples; especially sodium azide as this chemical would affect the enzymatic activity of the conjugate, generating false negative results.
3. Samples have to be clearly identified with codes or names in order to avoid misinterpretation of results. When the kit is used for the screening of blood units, bar code labelling and reading is strongly recommended.
4. Haemolysed (red) and visibly hyperlipemic ("milky") samples have to be discarded as they could generate false results. Samples containing residues of fibrin or heavy particles or microbial filaments and bodies should be discarded as they could give rise to false results.
5. Sera and plasma can be stored at +2...8°C for up to five days after collection. For longer storage periods, samples can be stored frozen at -20°C for several months. Any frozen samples should not be frozen/thawed more than once as this may generate particles that could affect the test result.
6. If particles are present, centrifuge at 2.000 rpm for 20 min or filter using 0.2-0.8µ filters to clean up the sample for testing.
7. Do not use heat inactivated samples as they could give origin to false reactivity.
8. Storage of diluted samples is not recommended and can adversely affect test performance.
9. Mix the samples before the use.

#### **H. PREPARATION OF COMPONENTS AND WARNINGS**

A study conducted on an opened kit has not pointed out any relevant loss of activity up to 2 months.

Allow the kit reagents to reach room temperature (+18...24°C) at least 30 minutes before use.

Take only the volume necessary for the testing. Return the unused portion at +2...8°C.

##### **1. Microplate**

Allow the microplate to reach room temperature (about 1 hr) before opening the container.

Unused strips have to be placed back into the aluminum pouch firmly zipped and stored at +2...8°C.

When opened the first time, residual strips are stable up to two months or until the humidity indicator inside the desiccant bag turns from yellow to green.

##### **2. Negative Control**

Ready to use. Mix well on vortex before use.

##### **3. Positive Controls**

Ready to use. Mix well on vortex before use. Handle this component as potentially infectious.

##### **4. Conjugate # 1**

Ready to use. Mix well on vortex before use.

Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.

If this component has to be transferred use only plastic, possibly sterile disposable containers.

## 5. Conjugate # 2 Concentrate 100x

The 100x concentrate reagent has to be diluted in Conjugate Diluent # 2.

Example: Add 250 µL of concentrated Conjugate # 2 to 24,75 mL of Conjugate # 2 Diluent.

Once diluted mix well end-over-end before use.

Once diluted, be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.

If this component has to be transferred use only plastic, possibly sterile disposable containers.

## 6. Conjugate # 2 Diluent

Ready to use. Mix well on vortex before use.

Reagent is used for the dilution of the Conjugate # 2

## 7. Sample Diluent

Ready to use reagent. Mix well on vortex before use.

## 8. Substrate TMB

Ready-to-use component.

Mix well on vortex before use.

Be careful not to contaminate the liquid with oxidizing chemicals, air-driven dust or microbes.

Do not expose to strong illumination, oxidizing agents and metallic surfaces.

If this component has to be transferred use only plastic, possible sterile disposable container.

## 9. Stop Solution

Ready to use. Mix well on vortex before use.

## 10. Wash Buffer Concentrate 25x

The 25x concentrated solution has to be diluted with EIA grade water and mixed gently end-over-end before use.

Example: Add 50 mL of 25x wash solution to 1200 mL of deionised water.

As some salt crystals may be present into the vial, take care to dissolve all the content during the preparation of the solution.

In the preparation avoid foaming as the presence of bubbles could give origin to a bad washing efficiency.

Once diluted, the wash solution is stable for 6 days at +2...8°C.

## I. INSTRUMENTS AND TOOLS USED IN COMBINATION WITH THE KIT

1. Micropipettes have to be calibrated to deliver the correct volume required by the assay and must be submitted to regular decontamination (household alcohol, 10% solution of bleach, hospital grade disinfectants) of those parts that could accidentally come in contact with the sample. They should also be regularly maintained. Decontamination of spills or residues of kit components should also be carried out regularly. They should also be regularly maintained in order to show a precision of 1% and a trueness of  $\pm 2\%$ .
2. The ELISA incubator has to be set at +37°C (tolerance of  $\pm 0.5^\circ\text{C}$ ) and regularly checked to ensure the correct temperature is maintained. Both dry incubators and water baths are suitable for the incubations, provided that the instrument is validated for the incubation of ELISA tests.
3. The ELISA washer is extremely important to the overall performances of the assay. The washer must

be carefully validated and correctly optimized using the kit controls and reference panels, before using the kit for routine laboratory tests.

**5-3-5 washing cycles (aspiration + dispensation of 400 µL/well of washing solution = 1 cycle) are sufficient to ensure that the assay performs as expected.** In order to set correctly their number, it is recommended to run an assay with the kit controls

and well characterized negative and positive reference samples, and check to match the values reported below in the section "Internal Quality control". Regular calibration of the volumes delivered by, and maintenance (decontamination and cleaning of needles) of the washer has to be carried out according to the instructions of the manufacturer.

4. Incubation times have a tolerance of  $\pm 5\%$  (or for 1<sup>st</sup> incubation tolerance between 57 min to 63 min, for 2<sup>nd</sup>, 3<sup>rd</sup> and 4<sup>nd</sup> incubation tolerance between 29 min to 31 min).
5. The ELISA microplate reader has to be equipped with a reading filter of 450nm and ideally with a second filter (600-650nm) for blanking purposes. Its standard performances should be (a) bandwidth  $\leq 10\text{nm}$ ; (b) absorbance range from 0 to  $\geq 2.0$ ; (c) linearity to  $\geq 2.0$ ; repeatability  $\geq 1\%$ . Blanking is carried out on the well identified in the section "Assay Procedure". The optical system of the reader has to be calibrated regularly to ensure that the correct optical density is measured. It should be regularly maintained according to the manufacturer 's instructions.

## L. PRE ASSAY CONTROLS AND OPERATIONS

1. Check the expiration date of the kit printed on the external label. Do not use the device if expired.
2. Check that the liquid components are not contaminated by visible particles or aggregates. Check that the Substrate is colourless or pale blue by aspirating a small volume of it with a sterile plastic pipette. Check that no breakage occurred in transportation and no spillage of liquid is present inside the box. Check that the aluminium pouch, containing the microplate, is not punctured or damaged.
3. Dilute the Conjugate # 2 as described above.
4. Dilute all the content of the 25x concentrated wash buffer as described above.
5. Allow all the other components to reach room temperature (about 1 hr) and then mix gently on vortex all liquid reagents as described.
6. Set the ELISA incubator at  $+37^\circ\text{C} \pm 0.1^\circ\text{C}$  and prepare the ELISA washer by priming with the diluted wash solution, according to the manufacturers instructions. Set the right number of washing cycles as found in the validation of the instrument for its use with the kit.
7. Check that the ELISA reader is turned on or ensure it will be turned on at least 20 minutes before reading.
8. If using an automated work station, turn on, check settings and be sure to use the right assay protocol.
9. Check that the micropipettes are set to the required volume.
10. Check that all the other equipment is available and ready to use.
11. In case of problems, do not proceed further with the test and advise the supervisor.

### M. ASSAY PROCEDURE (Manual)

The assay has to be carried out according to what reported below, taking care to maintain the same incubation time for all the samples in testing.

- Place the required number of microwells in the microwell holder. Leave the 1<sup>st</sup> well empty for the operation of blanking (optional).
- Dispense 100 µL of Sample Diluent to each well except the well for operation of blanking.
- Add 100 µL of Sample, 100 µL of Negative and Positive Controls in duplicate pipetting up and down for homogenization.  
Pipette gently avoiding overflowing and contaminating adjacent wells.  
Seal strips securely with microplate sealer.
- Incubate the microplate for **60 min at +37°C**.
- After incubation, remove solution from the wells by inverting the microplate and tapping dry on paper towel. Wash the microtiter plate 5 cycles according to the washing procedure (chapter I.3).
- Add 200 µL of Conjugate # 1 into each well, except the 1<sup>st</sup> blanking well, and cover with the sealer.  
**Important note:** Be careful not to touch the plastic inner surface of the well with the tip filled with the Conjugate. Contamination might occur.
- Incubate the microplate for **30 min at +37°C**.
- After incubation, remove solution from the wells by inverting the microplate and tapping dry on paper towel. Wash the microtiter plate 3 cycles according to the washing procedure (chapter I.3).
- Dispense 200 µL of diluted Conjugate # 2 into each well, except the 1<sup>st</sup> blanking well, and cover with the sealer.  
**Important note:** Be careful not to touch the plastic inner surface of the well with the tip filled with the Conjugate. Contamination might occur.
- Incubate the microplate for **30 min at +37°C**.
- After incubation, remove solution from the wells by inverting the microplate and tapping dry on paper towel. Wash the microtiter plate 5 cycles according to the washing procedure (chapter I.3).
- Pipette 150 µL of the Substrate TMB into each well, the blank well included.
- Incubate the microplate for **30 min at +18...25°C**.  
**Important note:** Do not expose to strong direct illumination. High background might be generated.
- Pipette 100 µL of Stop Solution into all the wells using the same pipetting sequence as in step 12 to stop the enzymatic reaction. Addition of acid will turn the positive control and positive samples from blue to yellow.
- Measure the color intensity of the solution in each well, as described in chapter I.5, at 450nm filter (reading) and possibly at 600-650nm (background subtraction), blanking the instrument on A1.  
Reading has to be carried out just after the addition of the Stop Solution and anyway not any longer than 30 minutes after its addition. Some self oxidation of the chromogen can occur leading to high background.

### Important notes:

- If the second filter is not available, ensure that no finger prints are present on the bottom of the microwell before reading at 450nm. Finger prints could generate false positive results on reading.
- Reading should ideally be performed immediately after the addition of the Stop Solution but definitely no longer than 20 minutes afterwards. Some self oxidation of the chromogen can occur leading to a higher background

### N. ASSAY SCHEME

STEP	PROCEDURE
Sample step	Add 100 µL of Sample Diluent. Add 100 µL of sample/control to the diluent and pipet up & down for mixing. Test Negative and Positive Controls in duplicate. Incubate 60 minutes at 37°C.
Wash step	Perform wash step 5x.
Conjugate 1	Add 200 µL of ready-to-use Conjugate 1 to wells of the microtiter plate. Incubate 30 minutes at 37°C.
Wash step	Perform wash step 3x.
Conjugate 2	Add 200 µL of diluted Conjugate 2 to wells of the microtiter plate. Incubate 30 minutes at 37°C
Wash step	Perform wash step 5x.
Substrate TMB	Add 150 µL of Substrate TMB to the wells of the microtiter plate.
Colour Development	Incubate 30 minutes at +18...25°C.
Stopping	Add 100 µL of Stop Solution and read OD at 450nm with reference wavelength at 600-650nm in an ELISA microplate reader.

An example of dispensation scheme is reported in the table below:

		Microplate											
		1	2	3	4	5	6	7	8	9	10	11	12
A	BLK	PC3	S8										
B	NC	S1	S9										
C	NC	S2											
D	PC1	S3											
E	PC1	S4											
F	PC2	S5											
G	PC2	S6											
H	PC3	S7											

Legenda: BLK = Blank NC = Negative Control  
PC1 = Positive Control HIV-1  
PC2 = Positive Control HIV-2  
PC3 = Positive Control HIV p-24  
S = Sample

### O. INTERNAL QUALITY CONTROL

A check is performed on the controls and or calibrator any time the kit is used in order to verify whether the expected OD450nm values have been matched in the analysis.

Ensure that the following parameters are met:

Check	Requirements
Blank well	< 0.100 OD450nm value
Negative Control (NC)	< 0.100 mean OD450nm value after blanking
Positive Control HIV-1	> 0.500 OD450nm value
Positive Control HIV-2	> 0.500 OD450nm value
Positive Control HIV p-24	> 0.500 OD450nm value

If the results of the test match the requirements stated above, proceed to the next section.

If they do not, do not proceed any further and perform the following checks:

Problem	Check
Blank well > 0.100 OD450nm	1. that the Substrate solution has not become contaminated during the assay
Negative Control (NC) > 0.100 OD450nm after blanking	1. that the washing procedure and the washer settings are as validated in the pre qualification study; 2. that the proper washing solution has been used and the washer has been primed with it before use; 3. that no mistake has been done in the assay procedure (dispensation of positive control instead of negative control); 4. that no contamination of the negative control or of the wells where the control was dispensed has occurred due to positive samples, to spills or to the enzyme conjugate; 5. that micropipettes have not become contaminated with positive samples or with the enzyme conjugate 6. that the washer needles are not blocked or partially obstructed.
Positive Controls < 0.500 OD450nm	1. that the procedure has been correctly performed; 2. that no mistake has occurred during the distribution of the control (dispensation of negative control instead of positive control). 3. that the washing procedure and the washer settings are as validated in the pre qualification study; 4. that no external contamination of the positive control has occurred.

If any of the above problems have occurred, report the problem to the supervisor for further actions.

## P. CALCULATION OF RESULTS

### P.1. Validity

Two Negative Controls (NC), two HIV-1 Positive Controls (PC1), two HIV-2 Positive Controls (PC2) and two HIV-1 p24 Controls (PC3) should be included in each run. The results for the controls should be within the acceptance criteria before any sample results can be interpreted.

Calculation of the Negative Control Mean  $NC_{MEAN}$ :

Example:

Absorbance NC  
0.021  
0.025  
-----  
0.046

$$NC_{MEAN} = 0.046 / 2 = 0.023$$

The mean of the absorbance of the Negative Controls, after blanking, must be less than 0.100. If the mean value is greater than or equal to 0.100, the run should be repeated.

Calculation of the HIV-1 Positive Control Mean  $HIV1-PC_{MEAN}$ :

Example:

Absorbance HIV1-PC  
1.545  
1.239  
-----  
2.784

$$HIV1-PC_{MEAN} = 2.784 / 2 = 1.392$$

The mean of the absorbance of the HIV-1 Positive Controls must be greater than 0.500. If the mean value is less than or equal to 0.500, the run should be repeated.

Calculation of the HIV-2 Positive Control Mean  $HIV2-PC_{MEAN}$ :

Example:

Absorbance HIV2-PC  
1.459  
1.343  
-----  
2.802

$$HIV2-PC_{MEAN} = 2.802 / 2 = 1.401$$

The mean of the absorbance of the HIV-2 Positive Controls must be greater than 0.500. If the mean value is less than or equal to 0.500, the run should be repeated.

Calculation of the HIV-1 p24 Positive Control Mean  $HIV-1\ p24-PC_{MEAN}$ :

Example:

Absorbance HIV-1 p24-PC  
2.223  
2.172  
-----  
4.395

$$HIV-1\ p24-PC_{MEAN} = 4.395 / 2 = 2.198$$

The mean of the absorbance of the HIV-1 p24 Positive Controls must be greater than 0.500. If the mean value is less than or equal to 0.500, the run should be repeated.

**P.2. Cut-Off Calculation**

The tests results are calculated by means of a cut-off value determined with the following formula on the mean OD<sub>450nm</sub> value of the Negative Control (NC) - OD<sub>450nm</sub> Blank:

$$(NC_{mean} - blank) + 0.170 = \text{Cut-Off (Co)}$$

The value found for the test is used for the interpretation of results as described in the next paragraph.

**Important note:** When the calculation of results is performed by the operating system of an ELISA automated work station, ensure that the proper formulation is used to calculate the cut-off value and generate the correct interpretation of results.

**Q. INTERPRETATION OF RESULTS**

1. Specimens with absorbance values less than the cut-off value are considered not reactive by the criteria of this immunoassay, and may be considered negative for antibodies to HIV1+2 and HIV-1 p24. Further testing is not required.
2. Specimens with absorbance values equal to or greater than the cut-off are considered to be reactive or positive for HIV-1 and/or HIV-2 antibodies or HIV-1 p24. These specimens (using the original specimen) should be re-tested in duplicate before final confirmation of the result.
3. Initially reactive specimens, which do not react in either of the duplicate repeat tests, are considered negative for antibodies of HIV1+2 and HIV-1 p24. Further testing is not required.
4. If one of both retested values is equal to or greater than the cut-off value, the specimen is considered repeatedly reactive. Specimens that have been found repeatedly reactive are interpreted to be positive for the presence of antibodies to HIV-1 and/or HIV-2 or HIV-1 p24. In most settings it is appropriate to investigate repeatedly reactive specimens by additional, more specific tests.

**Important notes:**

1. Interpretation of results should be done under the supervision of the responsible of the laboratory to reduce the risk of judgment errors and misinterpretations.
2. Repeatedly reactive specimens should be submitted to a Confirmation Assay before diagnosis of HIV infection is released.
3. When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.
4. Diagnosis of HIV infection has to be done and released to the patient only by a qualified medical doctor.

**R. PERFORMANCE CHARACTERISTICS**

Evaluation of Performances has been conducted in accordance to what reported in the Common Technical Specifications (CTS) as required by art. 5, Chapter 3 of IVD Directive 98/79/EC.

The performances evaluation was carried out both in an external centre and in Adaltis' laboratories as well to complete the study.

**R.1. SPECIFICITY AND SENSITIVITY**

**R.1.1 Specificity**

The specificity evaluated on:

- ❖ 5014 European unselected blood donors from two centres including 40 samples from first time donors
- ❖ 201 Hospitalized patients
- ❖ 116 Potentially interfering samples:
  - ✓ 27 Pregnancy samples (6 monopar, 17 multipar and 4 unknown)
  - ✓ 27 Rheumatoid Factor (RF) positive samples
  - ✓ 25 Hemolyzed samples including heavy hemolyzed
  - ✓ 27 with other infection (4 HBV, 5 EBV, 5 HCV, 4 HSV, 5 Rubella, 4 Toxoplasmosis)
  - ✓ 5 Hyperbilirubinemia samples
  - ✓ 5 Hyperlipemia samples
- ❖ 50 Samples [25 Positives (HIV-1, HIV-2 and HIV-1 p24) and 25 Negative] have been tested and no difference due to the method of sample preparation (plasma-serum, citrate, EDTA and Heparin) has been observed.

From the 5014 European unselected blood donors, 19 (nineteen) samples were initially reactive (IR) and tested in duplicate. These resulted in 17 (seventeen) repeated reactive (RR) and they were tested with a confirmatory test. One sample remained indeterminate.

Actual Status	Test Results			
	IR	RR	NEG	Total
Negative	18	16	4997	5013
Positive	0	0	0	0
Indetermined	1	1	0	1
Total	19	17	4997	5014

A final specificity has been found at 99.7%.

**R.1.2 Analytical Sensitivity (Limit of Detection)**

The analytical sensitivity was determined by means of the WHO international standard for HIV-1 p24 Antigen, First International Reference Reagent, NIBSC code 90/636.

A dilution into a negative matrix was prepared (0.1-10 IU/mL). This dilution was tested in three different test-runs (three microplates) over two tests (two days). For each individual run, a new dilution of samples was made. The summary results are shown in Table below.

Run	1	2	3
HIV-1 p24 (IU/mL)	OD/CO		
10	11.9	12.6	12.3
5	8.0	10.7	9.8
2.5	4.3	6.2	5.7
1	1.8	2.8	2.6
0.5	1.0	1.6	1.4
0.25	0.6	1.0	0.8
0.1	0.4	0.5	0.5
0	0.2	0.3	0.3

The sensitivity limit has been estimated < 1.0 IU/mL

**R.1.3 Diagnostic Sensitivity**

The diagnostic sensitivity of the EIAGEN Detect HIV 4 Total Screening Kit was based on testing of a panel of 500 HIV1/2 positive samples, 31 seroconversion panels

(supplied by BBI-Seracare/ZeptoMetrix) and 40 early seroconversion samples.

In details:

- ❖ 500 Positive Samples (400 HIV-1 and 100 HIV-2 including subtype A, B, C, D, F, G, H, J, K, CRF);
- ❖ 50 cell culture supernatants including different HIV-1 subtypes (group M subtype A to K, CRF01\_AE, CRF02\_AG, group N, group O) and HIV-2;
- ❖ 50 HIV-1 p24 Antigen Positive Samples
- ❖ 31 Seroconversion Panels (see table below)

Seroconversion Panel	Adaltis EIAgen	4 <sup>th</sup> gen. assay	3 <sup>rd</sup> gen. assay	HIV -1 p24 antigen	NAT
	First sample detected positive in the panel				
Y-PRB925	5		5	5	5
AB/PRB927	2	2	2	2	2
AC/PRB928	2	2	2	2	2
AD/PRB929	3	3	6	3	3
AE PRB930	1		3	1	1
PRB933	2	2	2	2	2
PRB934	1	1	1	1	1
AP/PRB940	1	2	3	2	1
AQ/PRB941	3	4	4	3	2
AS/PRB943	3	3	6	3	2
AU PRB945	3		4	3	1
AW/PRB947	2	2	2	2	1
AX PRB948	3	4	4+	3	3
AZPRB950	2		4	2	2
BA PRB 951	3	3	5	3	3
PRB952	3	3	4	3	2
BE/PRB955	2	2	4	2	1
BF PRB 956	4	4	5	4	2
BI/PRB959	1	1	3	1	1
PRB960	7	8	9+	8	8
PRB962	5	5	6+	5	3
PRB963	6	6	7	6	5
PRB964	6+	6	ND	6	4
PRB965	4	2	4	2	1
PRB967	4	4	4	4	2
PRB968	7	7	7	7	5
6240	7	8	9	8	6
9028	6	6	ND	6	6
9077	12	12	14	12	12
9079	9	9	11	9	8
12007	4	4	5	ND	4

- ❖ 42 Early Seroconversion Samples (positive for HIV-1 p24 core antigen and HIV antibodies are absent or weakly present - indeterminate result on Western Blot test).

A final sensitivity has been found 100%.

## R.2. PRECISION

The precision of the device was assessed by determining its values in a within and between runs. In the tables below results are reported for a Negative sample and Positive samples.

### R.2.1 Intra lot results:

*EIAgen Detect HIV 4 Total Screening code: 081311*

1<sup>st</sup> Lot

Sample	S/Co Mean	Precision - %CV		
		Within run	Between Run	Total
Negative	0.37	23.2	4.5	23.6
Pos A-HIV 1	1.80	10.6	9.4	14.1
Pos A-HIV 2	1.91	10.7	7.7	13.1
Pos HIV-1 p24	2.14	10.9	4.6	11.8

*EIAgen Detect HIV 4 Total Screening code: 081311*

2<sup>nd</sup> Lot

Sample	S/Co Mean	Precision - %CV		
		Within run	Between Run	Total
Negative	0.40	12.6	10.8	16.6
Pos A-HIV 1	1.47	6.9	7.4	10.3
Pos A-HIV 2	1.53	7.6	8.1	11.1
Pos HIV-1 p24	1.78	4.6	4.9	6.7

*EIAgen Detect HIV 4 Total Screening code: 081311*

3<sup>th</sup> Lot

Sample	S/Co Mean	Precision - %CV		
		Within run	Between Run	Total
Negative	0.52	18.7	5.4	19.5
Pos A-HIV 1	1.65	8.6	4.6	9.8
Pos A-HIV 2	1.78	9.7	4.9	10.9
Pos HIV-1 p24	2.05	7.4	4.6	8.7

*EIAgen Detect HIV 4 Total Screening code: 081311*

4<sup>th</sup> Lot

Sample	S/Co Mean	Precision - %CV		
		Within run	Between Run	Total
Negative	0.58	19.6	3.0	19.8
Pos A-HIV 1	1.59	9.4	1.4	9.5
Pos A-HIV 2	1.66	11.4	2.64	11.7
Pos HIV-1 p24	2.17	9.1	2.0	9.3

### R.2.2 Inter lot results:

Sample	S/Co Mean	Precision - %CV	
		Between run	Total
Negative	0.48	18.4	26.5
Pos A-HIV 1	1.63	9.7	13.3
Pos A-HIV 2	1.72	9.9	14.2
Pos HIV-1 p24	2.03	8.6	12.1

### R.2.3 Specifications:

Intra Lot:

Within Run:

%CV on positive samples < 15%

%CV on negative sample < 25%

Between Run:

%CV on positive samples < 15%

%CV on negative sample < 25%

Inter lot:

Between Lot:

%CV on positive samples < 15%

%CV on negative sample < 25%

Total precision:

%CV on positive samples < 20%

%CV on negative sample < 30%

## S. SUGGESTIONS FOR TROUBLESHOOTING

Adherence to assay procedure and specifications, as well as a correct use of reagents and proper pipetting, may help to avoid the following kinds of errors.

ERROR	POSSIBLE CAUSES / SUGGESTIONS
OD very different ( $\pm 50\%$ ) from OD reported on QC	<ul style="list-style-type: none"> <li>- incorrect dispensing volume of reagents (suggestion: check the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes)</li> <li>-incorrect temperature or incorrect incubation time (suggestion: more care in the incubator maintenance; note down the beginning of the incubation)</li> <li>-error in washing or in photometer reading (suggestion: check operating or settings of respective instruments)</li> <li>-contamination of Substrate or Conjugate (suggestion: use only disposable and clean plastic containers)</li> </ul>
Low reproducible results	<ul style="list-style-type: none"> <li>-not constant dispensing volume of samples or reagents (suggestion: check the pipettes precision and the correspondence between the volume dispensed by the pipette and the one required by the assay; re-calibrate again pipettes)</li> <li>-error in washing or in reading (suggestion: check operating or settings of respective instruments)</li> <li>-contamination of Substrate (suggestion: use only disposable and clean plastic containers)</li> <li>-pollution or degradation of reagents (suggestion: use appropriate tips, disposable and clean plastic containers for reagents and high quality distilled or equivalent water)</li> </ul>
no colorimetric reaction after addition of substrate	<ul style="list-style-type: none"> <li>-some reagent not pipetted</li> <li>- strong contamination of Conjugates or Substrate</li> <li>-errors in performing the assay procedure (e.g. accidental pipetting of reagents in a wrong sequence or from the wrong vial, etc.)</li> </ul>
too low reaction (too low ODs)	<ul style="list-style-type: none"> <li>-incubation time too short, incubation temperature too low</li> <li>-incorrect conjugate dilution</li> </ul>
too high reaction (too high ODs)	<ul style="list-style-type: none"> <li>-incorrect conjugate dilution</li> <li>-incubation time too long, incubation temperature too high</li> <li>-water quality for wash buffer insufficient (low grade of deionization)</li> <li>-insufficient washing (conjugates not properly removed)</li> </ul>
unexplainable outliers	<ul style="list-style-type: none"> <li>-contamination of pipettes, tips or containers</li> <li>-inconstant and insufficient washing (conjugates not properly removed)</li> </ul>
too high within-run CV%	<ul style="list-style-type: none"> <li>-reagents and/or strips not pre-warmed to Room Temperature prior to use</li> <li>- plate washer is not washing correctly (suggestion: clean washer head)</li> </ul>

ERROR	POSSIBLE CAUSES / SUGGESTIONS
too high between-run CV%	<ul style="list-style-type: none"> <li>-incubation conditions not constant (time, temperature)</li> <li>-controls and samples not dispensed at the same time (with the same intervals) (check pipetting order)</li> <li>-person-related variation</li> </ul>

## T. AUTOMATION

The procedure identified in this instruction for use is for manual testing only. When using automated instruments, follow the procedures that are contained in the operator's manual provided by the device manufacturer. Laboratories must follow their approved validation procedures to demonstrate compatibility of this product on automated systems.

## U. LIMITATIONS

- The user of this kit is advised to carefully read and understand the instructions for use. Strict adherence to the protocol is necessary to obtain reliable test results. In particular, correct sample and reagent pipetting, along with careful washing and timing of incubation steps is essential for accurate, reproducible detection of HIV-1 and HIV-2 antibodies and p24 antigens.
- If possible, use fresh serum or plasma samples. Sample degradation as well as multiple freeze-thaw cycles may cause erroneous results. **Do not use heat-inactivated samples.**
- Falsely reactive test results can be expected with a test kit of this nature. The proportion of reactivities will depend on the sensitivity and specificity of the test kit and on the prevalence of HIV-1 and HIV-2 antibodies in the population to be screened.
- After the EIAgen Detect HIV 4 Total Screening Kit is performed, repeatedly reactive samples should be submitted for additional testing using Western Blot (WB), Indirect Immunofluorescence Assay (IFA) or Radioimmunoassay (RIA) tests. The determination that a person's samples contains antibodies to HIV and/or HIV antigens has extensive medical, social, psychological and economic implications. It is recommended that confidentiality, appropriate counselling and medical evaluation be considered an essential aspect of the testing sequence.
- AIDS and AIDS-related conditions are clinical diseases and their diagnosis can only be established clinically. EIA testing alone cannot be used to diagnose AIDS. A non-reactive test result at any point in the testing sequence does not preclude the possibility of exposure to or infection with HIV. The risk of an asymptomatic person, who is repeatedly reactive, of developing AIDS and/or AIDS-related conditions, is not known.
- The test result should be used in conjunction with all other clinical and diagnostic data.
- Antibodies to HIV may occur due to voluntary participation in an HIV vaccine study. Interpretation of this diagnostic test will depend on the type of vaccine given. Correlation with the medical history and additional testing may be necessary to accurately diagnose HIV in vaccine volunteers.

## BIBLIOGRAPHY

1. Alizon, M., Sonigo, P., Barré-Sinoussi, F., Chermann, J.-C., Tiollais, P., Montagnier, L. and Wain-Hobson, S., 1984. Molecular Cloning of Lymphadenopathy-Associated Virus. *Nature* 312:757-760.
2. Barré-Sinoussi, F., Chermann, J.-C., Rey, F., Nugeyre, M.T., Chamaret, S., Gruest, J., Dauguet, C., Axler-Blin, C., Vézinet-Brun, F., Rouzioux, C., Rozenbrum, W. and Montagnier, L. 1983. Isolation of a T-lymphotropic Retrovirus from a Patient at Risk for Acquired Immune Deficiency Syndrome (AIDS). *Science* 220:868-871.
3. Clavel, F., Mansinho, K., Chamaret, S. et al. 1987. Human Immunodeficiency Virus Type 2 Infection Associated with AIDS in West Africa. *N. Engl. J. Med.* 316:1180-1185.
4. Gallo, R.C., Salahuddin, S.Z., Popovic, M., Shearer, G.M., Kaplan, M., Haynes, B.F., Palker, T.J., Redfield, R., Oleske, J., Safal, B., White, G., Foster, P. and Markham, P.D. 1984. Frequent Detection and Isolation of Cytopathic Retroviruses (HTLV-III) from Patients with AIDS and at Risk for AIDS. *Science* 224:500-503.
5. Gold, J. and Dwyer, J., 1994. A Short History of AIDS. *Med. J. Aust.* 160:251-252.
6. Hahn, B.N., Shaw, G.M., Arya, S.K., Popovic, M., Gallo, R.C. and Wong-Staal, F., 1984. Molecular Cloning and Characterization of the HTLV-III Virus Associated with AIDS. *Nature* 312:166-169.
7. IVD Directive 98/79/CE, Common Technical Specifications (CTS) – Annex II, List A.
8. Lin, H.J. 1995. Laboratory Tests for Human Immunodeficiency Viruses. *J. Int. Fed. Clin. Chem.* 7:61-65.
9. Luciw, P.A., Potter, S.J., Steimer, K., Dina, D. and Levy, J.A., 1984. Molecular Cloning of AIDS-Associated Retrovirus. *Nature* 312:760-763.
10. Ly, T.D., Laperche, S., Brennan, C., Vallari, A., Ebel, A., Hunt, J., Martin, L., Daghfal, D., Schochetman, G. And Devare, S. 2004. Evaluation of the sensitivity and specificity of six HIV combined p24 antigen and antibody assays. *J. Virol. Meth.* 122: 185-194.
11. Popovic, M., Sarngadharan, M.G, Read, E., and Gallo, R.C., 1984. Detection, Isolation, and Continuous Production of Cytopathic Retrovirus (HTLV-III) from Patient with AIDS and Pre-AIDS. *Science* 224:497-500.
12. Sarngadharan, M.G., Popovic, M., Bruch, L., Schüpbach, J. and Gallo, R.C., 1984. Antibodies Reactive with Human T-Lymphotropic Retroviruses (HTLV-III) in the Serum of Patients with AIDS. *Science* 224:506-508.
13. Saville, R.D., Constantine, N.T., Cleghorn, F.R., Jack, N., Bartholomew, C., Edwards, J., Gomez, P. and Blattner, W.A. 2001. Fourth-generation enzyme-linked immunosorbent assay for the simultaneous detection of human immunodeficiency virus antigen and antibody. *J. Clin. Microbiol.* 39 (7): 2518-2524.
14. Sehulster, L.M., Hollinger, F.B., Dreesman, G.R. and Melnick, J.L., 1981. Immunological and Biophysical Alteration of Hepatitis B Virus Antigens by Sodium Hypochlorite Disinfection, *App. and Environ. Microbiol.* 42:762-767.
15. Sinicco, A., For a, R., Scalandra, M., Lucchini, A., Caramello, P. and Giovanni, P. 1993. Risk of Developing AIDS after Primary Acute HIV-1 Infection. *J. Acquir. Immune. Defic. Syndr.* 6:575-581.
16. Spire, B., Montagnier, L., Barré-Sinoussi, F. and Chermann, J.-C., 1984. Inactivation of Lymphadenopathy Associated Virus by Chemical Disinfectants. *Lancet*: 889-901, Oct. 20.
17. Vézinet-Brun, F., Barré-Sinoussi, F., Salmot, A.G., Christol, D. Montagnier, L., Rouzioux, C., Klatzmann, D., Rozenbaum, W., Gluckmann, J.C. and Chermann, J.-C., 1984. Detection of IgG Antibodies to Lymphadenopathy-Associated Virus in Patients with AIDS or Lymphadenopathy Syndrome. *Lancet*: 1253-1256, June 9.
18. Weber, B., Thorstensson, R., Tanprasert, S., Schmitt, U. And Melchior, W. 2003. Reduction of the diagnostic window in three cases of human immunodeficiency-1 subtype E primary infection with fourth-generation HIV screening assays. *Vox Sanguinis* 85: 73-79.
19. World Health Organization. 2004. HIV assays: operational characteristics (Phase 1): report 15 antigen/antibody ELISAs. [www.who.int/diagnostics\\_laboratory/publications/en/HIV\\_Report15.pdf](http://www.who.int/diagnostics_laboratory/publications/en/HIV_Report15.pdf).

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
## BIBLIOGRAFIA

1. Alizon, M., Sonigo, P., Barré-Sinoussi, F., Chermann, J.-C., Tiollais, P., Montagnier, L. and Wain-Hobson, S., 1984. Molecular Cloning of Lymphadenopathy-Associated Virus. *Nature* 312:757-760.
2. Barré-Sinoussi, F., Chermann, J.-C., Rey, F., Nugeyre, M.T., Chamaret, S., Gruest, J., Daugey, C., Axler-Blin, C., Vézinet-Brun, F., Rouzioux, C., Rozenbrum, W. and Montagnier, L. 1983. Isolation of a T-lymphotropic Retrovirus from a Patient at Risk for Acquired Immune Deficiency Syndrome (AIDS). *Science* 220:868-871.
3. Clavel, F., Mansinho, K., Chamaret, S. et al. 1987. Human Immunodeficiency Virus Type 2 Infection Associated with AIDS in West Africa. *N. Engl. J. Med.* 316:1180-1185.
4. Gallo, R.C., Salahuddin, S.Z., Popovic, M., Shearer, G.M., Kaplan, M., Haynes, B.F., Palker, T.J., Redfield, R., Oleske, J., Safal, B., White, G., Foster, P. and Markham, P.D. 1984. Frequent Detection and Isolation of Cytopathic Retroviruses (HTLV-III) from Patients with AIDS and at Risk for AIDS. *Science* 224:500-503.
5. Gold, J. and Dwyer, J., 1994. A Short History of AIDS. *Med. J. Aust.* 160:251-252.
6. Hahn, B.N., Shaw, G.M., Arya, S.K., Popovic, M., Gallo, R.C. and Wong-Staal, F., 1984. Molecular Cloning and Characterization of the HTLV-III Virus Associated with AIDS. *Nature* 312:166-169.
7. IVD Directive 98/79/CE, Common Technical Specifications (CTS) – Annex II, List A.
8. Lin, H.J. 1995. Laboratory Tests for Human Immunodeficiency Viruses. *J. Int. Fed. Clin. Chem.* 7:61-65.
9. Luciw, P.A., Potter, S.J., Steimer, K., Dina, D. and Levy, J.A., 1984. Molecular Cloning of AIDS-Associated Retrovirus. *Nature* 312:760-763.
10. Ly, T.D., Laperche, S., Brennan, C., Vallari, A., Ebel, A., Hunt, J., Martin, L., Daghfal, D., Schochetman, G. And Devare, S. 2004. Evaluation of the sensitivity and specificity of six HIV combined p24 antigen and antibody assays. *J. Virol. Meth.* 122: 185-194.
11. Popovic, M., Sarngadharan, M.G., Read, E., and Gallo, R.C., 1984. Detection, Isolation, and Continuous Production of Cytopathic Retrovirus (HTLV-III) from Patient with AIDS and Pre-AIDS. *Science* 224:497-500.
12. Sarngadharan, M.G., Popovic, M., Bruch, L., Schüpbach, J. and Gallo, R.C., 1984. Antibodies Reactive with Human T-Lymphotropic Retroviruses (HTLV-III) in the Serum of Patients with AIDS. *Science* 224:506-508.
13. Saville, R.D., Constantine, N.T., Cleghorn, F.R., Jack, N., Bartholomew, C., Edwards, J., Gomez, P. and Blattner, W.A. 2001. Fourth-generation enzyme-linked immunosorbent assay for the simultaneous detection of human immunodeficiency virus antigen and antibody. *J. Clin. Microbiol.* 39 (7): 2518-2524.
14. Sehulster, L.M., Hollinger, F.B., Dreesman, G.R. and Melnick, J.L., 1981. Immunological and Biophysical Alteration of Hepatitis B Virus Antigens by Sodium Hypochlorite Disinfection, *App. and Environ. Microbiol.* 42:762-767.
15. Sinicco, A., For a, R., Scalandra, M., Lucchini, A., Caramello, P. and Giovanni, P. 1993. Risk of Developing AIDS after Primary Acute HIV-1 Infection. *J. Acquir. Immune. Defic. Syndr.* 6:575-581.
16. Spire, B., Montagnier, L., Barré-Sinoussi, F. and Chermann, J.-C., 1984. Inactivation of Lymphadenopathy Associated Virus by Chemical Disinfectants. *Lancet*: 889-901, Oct. 20.
17. Vézinet-Brun, F., Barré-Sinoussi, F., Salmot, A.G., Christol, D. Montagnier, L., Rouzioux, C., Klatzmann, D., Rozenbaum, W., Gluckmann, J.C. and Chermann, J.-C., 1984. Detection of IgG Antibodies to Lymphadenopathy-Associated Virus in Patients with AIDS or Lymphadenopathy Syndrome. *Lancet*: 1253-1256, June 9.
18. Weber, B., Thorstensson, R., Tanprasert, S., Schmitt, U. And Melchior, W. 2003. Reduction of the diagnostic window in three cases of human immunodeficiency-1 subtype E primary infection with fourth-generation HIV screening assays. *Vox Sanguinis* 85: 73-79.
19. World Health Organization. 2004. HIV assays: operational characteristics (Phase 1): report 15 antigen/antibody ELISAs. [www.who.int/diagnostics\\_laboratory/publications/en/HIV\\_Report15.pdf](http://www.who.int/diagnostics_laboratory/publications/en/HIV_Report15.pdf).

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## IRODALOMJEGYZÉK

1. Alizon, M., Sonigo, P., Barré-Sinoussi, F., Chermann, J.-C., Tiollais, P., Montagnier, L. and Wain-Hobson, S., 1984. Molecular Cloning of Lymphadenopathy-Associated Virus. *Nature* 312:757-760.
2. Barré-Sinoussi, F., Chermann, J.-C., Rey, F., Nugeyre, M.T., Chamaret, S., Gruest, J., Dauguet, C., Axler-Blin, C., Vézinet-Brun, F., Rouzioux, C., Rozenbrum, W. and Montagnier, L. 1983. Isolation of a T-lymphotropic Retrovirus from a Patient at Risk for Acquired Immune Deficiency Syndrome (AIDS). *Science* 220:868-871.
3. Clavel, F., Mansinho, K., Chamaret, S. et al. 1987. Human Immunodeficiency Virus Type 2 Infection Associated with AIDS in West Africa. *N. Engl. J. Med.* 316:1180-1185.
4. Gallo, R.C., Salahuddin, S.Z., Popovic, M., Shearer, G.M., Kaplan, M., Haynes, B.F., Palker, T.J., Redfield, R., Oleske, J., Safal, B., White, G., Foster, P. and Markham, P.D. 1984. Frequent Detection and Isolation of Cytopathic Retroviruses (HTLV-III) from Patients with AIDS and at Risk for AIDS. *Science* 224:500-503.
5. Gold, J. and Dwyer, J., 1994. A Short History of AIDS. *Med. J. Aust.* 160:251-252.
6. Hahn, B.N., Shaw, G.M., Arya, S.K., Popovic, M., Gallo, R.C. and Wong-Staal, F., 1984. Molecular Cloning and Characterization of the HTLV-III Virus Associated with AIDS. *Nature* 312:166-169.
7. IVD Directive 98/79/CE, Common Technical Specifications (CTS) – Annex II, List A.
8. Lin, H.J. 1995. Laboratory Tests for Human Immunodeficiency Viruses. *J. Int. Fed. Clin. Chem.* 7:61-65.
9. Luciw, P.A., Potter, S.J., Steimer, K., Dina, D. and Levy, J.A., 1984. Molecular Cloning of AIDS-Associated Retrovirus. *Nature* 312:760-763.
10. Ly, T.D., Laperche, S., Brennan, C., Vallari, A., Ebel, A., Hunt, J., Martin, L., Daghfal, D., Schochetman, G. and Devare, S. 2004. Evaluation of the sensitivity and specificity of six HIV combined p24 antigen and antibody assays. *J. Virol. Meth.* 122: 185-194.
11. Popovic, M., Sarngadharan, M.G., Read, E., and Gallo, R.C., 1984. Detection, Isolation, and Continuous Production of Cytopathic Retrovirus (HTLV-III) from Patient with AIDS and Pre-AIDS. *Science* 224:497-500.
12. Sarngadharan, M.G., Popovic, M., Bruch, L., Schüpbach, J. and Gallo, R.C., 1984. Antibodies Reactive with Human T-Lymphotropic Retroviruses (HTLV-III) in the Serum of Patients with AIDS. *Science* 224:506-508.
13. Saville, R.D., Constantine, N.T., Cleghorn, F.R., Jack, N., Bartholomew, C., Edwards, J., Gomez, P. and Blattner, W.A. 2001. Fourth-generation enzyme-linked immunosorbent assay for the simultaneous detection of human immunodeficiency virus antigen and antibody. *J. Clin. Microbiol.* 39 (7): 2518-2524.
14. Sehulster, L.M., Hollinger, F.B., Dreesman, G.R. and Melnick, J.L., 1981. Immunological and Biophysical Alteration of Hepatitis B Virus Antigens

- by Sodium Hypochlorite Disinfection, *App. and Environ. Microbiol.* 42:762-767.
15. Sinicco, A., For a, R., Scalandra, M., Lucchini, A., Caramello, P. and Giovanni, P. 1993. Risk of Developing AIDS after Primary Accute HIV-1 Infection. *J. Acquir. Immune. Defic. Syndr.* 6:575-581.
  16. Spire, B., Montagnier, L., Barré-Sinoussi, F. and Chermann, J.-C., 1984. Inactivation of Lymphadenopathy Associated Virus by Chemical Dinsinfetants. *Lancet*: 889-901, Oct. 20.
  17. Vézinnet-Brun, F., Barré-Sinoussi, F., Salmot, A.G., Christol, D. Montagnier, L., Rouzioux, C., Klatzmann, D., Rozenbaum, W., Gluckmann, J.C. and Chermann, J.-C., 1984. Detection of IgG Antibodies to Lymphadenopathy-Associated Virus in Patients with AIDS or Lymphadenopathy Syndrome. *Lancet*: 1253-1256, June 9.
  18. Weber, B., Thorstensson, R., Tanprasert, S., Schmitt, U. And Melchior, W. 2003. Reduction of the diagnostic window in three cases of human immunodeficiency-1 subtype E primary infection with fourth-generation HIV screening assays. *Vox Sanguinis* 85: 73-79.
  19. World Health Organization. 2004. HIV assays: operational characteristics (Phase 1): report 15 antigen/antibody ELISAs. [www.who.int/diagnostics\\_laboratory/publications/en/HIV\\_Report15.pdf](http://www.who.int/diagnostics_laboratory/publications/en/HIV_Report15.pdf).

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
## BIBLIOGRAPHY

1. Alizon, M., Sonigo, P., Barré-Sinoussi, F., Chermann, J.-C., Tiollais, P., Montagnier, L. and Wain-Hobson, S., 1984. Molecular Cloning of Lymphadenopathy-Associated Virus. *Nature* 312:757-760.
2. Barré-Sinoussi, F., Chermann, J.-C., Rey, F., Nugeyre, M.T., Chamaret, S., Gruest, J., Daugey, C., Axler-Blin, C., Vézinet-Brun, F., Rouzioux, C., Rozenbrum, W. and Montagnier, L. 1983. Isolation of a T-lymphotropic Retrovirus from a Patient at Risk for Acquired Immune Deficiency Syndrome (AIDS). *Science* 220:868-871.
3. Clavel, F., Mansinho, K., Chamaret, S. et al. 1987. Human Immunodeficiency Virus Type 2 Infection Associated with AIDS in West Africa. *N. Engl. J. Med.* 316:1180-1185.
4. Gallo, R.C., Salahuddin, S.Z., Popovic, M., Shearer, G.M., Kaplan, M., Haynes, B.F., Palker, T.J., Redfield, R., Oleske, J., Safal, B., White, G., Foster, P. and Markham, P.D. 1984. Frequent Detection and Isolation of Cytopathic Retroviruses (HTLV-III) from Patients with AIDS and at Risk for AIDS. *Science* 224:500-503.
5. Gold, J. and Dwyer, J., 1994. A Short History of AIDS. *Med. J. Aust.* 160:251-252.
6. Hahn, B.N., Shaw, G.M., Arya, S.K., Popovic, M., Gallo, R.C. and Wong-Staal, F., 1984. Molecular Cloning and Characterization of the HTLV-III Virus Associated with AIDS. *Nature* 312:166-169.
7. IVD Directive 98/79/CE, Common Technical Specifications (CTS) – Annex II, List A.
8. Lin, H.J. 1995. Laboratory Tests for Human Immunodeficiency Viruses. *J. Int. Fed. Clin. Chem.* 7:61-65.
9. Luciw, P.A., Potter, S.J., Steimer, K., Dina, D. and Levy, J.A., 1984. Molecular Cloning of AIDS-Associated Retrovirus. *Nature* 312:760-763.
10. Ly, T.D., Laperche, S., Brennan, C., Vallari, A., Ebel, A., Hunt, J., Martin, L., Daghfal, D., Schochetman, G. And Devare, S. 2004. Evaluation of the sensitivity and specificity of six HIV combined p24 antigen and antibody assays. *J. Virol. Meth.* 122: 185-194.
11. Popovic, M., Sarngadharan, M.G, Read, E., and Gallo, R.C., 1984. Detection, Isolation, and Continuous Production of Cytopathic Retrovirus (HTLV-III) from Patient with AIDS and Pre-AIDS. *Science* 224:497-500.
12. Sarngadharan, M.G., Popovic, M., Bruch, L., Schüpbach, J. and Gallo, R.C., 1984. Antibodies Reactive with Human T-Lymphotropic Retroviruses (HTLV-III) in the Serum of Patients with AIDS. *Science* 224:506-508.
13. Saville, R.D., Constantine, N.T., Cleghorn, F.R., Jack, N., Bartholomew, C., Edwards, J., Gomez, P. and Blattner, W.A. 2001. Fourth-generation enzyme-linked immunosorbent assay for the simultaneous detection of human immunodeficiency virus antigen and antibody. *J. Clin. Microbiol.* 39 (7): 2518-2524.
14. Sehulster, L.M., Hollinger, F.B., Dreesman, G.R. and Melnick, J.L., 1981. Immunological and Biophysical Alteration of Hepatitis B Virus Antigens by Sodium Hypochlorite Disinfection, *App. and Environ. Microbiol.* 42:762-767.
15. Sinicco, A., For a, R., Scalandra, M., Lucchini, A., Caramello, P. and Giovanni, P. 1993. Risk of Developing AIDS after Primary Acute HIV-1 Infection. *J. Acquir. Immune. Defic. Syndr.* 6:575-581.
16. Spire, B., Montagnier, L., Barré-Sinoussi, F. and Chermann, J.-C., 1984. Inactivation of Lymphadenopathy Associated Virus by Chemical Disinfectants. *Lancet*: 889-901, Oct. 20.
17. Vézinet-Brun, F., Barré-Sinoussi, F., Salmot, A.G., Christol, D. Montagnier, L., Rouzioux, C., Klatzmann, D., Rozenbaum, W., Gluckmann, J.C. and Chermann, J.-C., 1984. Detection of IgG Antibodies to Lymphadenopathy-Associated Virus in Patients with AIDS or Lymphadenopathy Syndrome. *Lancet*: 1253-1256, June 9.
18. Weber, B., Thorstensson, R., Tanprasert, S., Schmitt, U. And Melchior, W. 2003. Reduction of the diagnostic window in three cases of human immunodeficiency-1 subtype E primary infection with fourth-generation HIV screening assays. *Vox Sanguinis* 85: 73-79.
19. World Health Organization. 2004. HIV assays: operational characteristics (Phase 1): report 15 antigen/antibody ELISAs. [www.who.int/diagnostics\\_laboratory/publications/en/HIV\\_Report15.pdf](http://www.who.int/diagnostics_laboratory/publications/en/HIV_Report15.pdf).